

CARE & CLEANING GUIDE

INSTRUCTIONS FOR USE

[PHILLIPS-SAFETY.COM](https://phillips-safety.com)

COPYRIGHT ©2026. ALL RIGHTS RESERVED



PHILLIPS SAFETY
PRODUCTS INC.

Your OrthoVest is built for daily use in the procedure suite. A few minutes of routine care keeps it clean, comfortable, and working properly.

DAILY (AFTER EACH USE)

- Remove your lead apron before taking off the OrthoVest.
- Wipe the frame and belt with a hospital-grade disinfectant wipe (quaternary ammonium or 70% isopropyl alcohol).
- Allow all parts to air dry before storing or putting back on.

CLEANING THE FRAME

- The frame is high-density polyethylene (HDPE) and resists most hospital disinfectants.
- Wipe all surfaces, including inside the keyhole slots.
- For heavier soil, wash with mild soap and warm water, rinse, and dry with a soft cloth.
- Do not use abrasive pads or scrapers, which can scratch the surface.

CLEANING THE PADDED BELT

- Wipe the fabric with a disinfectant wipe or a damp cloth with mild soap.
- Spot clean as needed; do not soak the padding.
- Air dry fully. Do not machine wash or machine dry.
- Avoid bleach on the fabric and webbing, as it can weaken and discolor the material.

HEAT PRECAUTIONS

- Do not autoclave, steam sterilize, or machine dry any part of the OrthoVest.
- Keep the frame away from radiators, heaters, and hot car interiors. High heat can warp the frame and change its fit.

WEEKLY INSPECTION

- Frame: check for cracks, stress marks, or warping, especially around keyhole slots.
- Studs: confirm each stud is tight and seats firmly in the keyhole slots.
- Buckles and adjusters: confirm they latch and release smoothly.
- Webbing: check for fraying, cuts, or loose stitching.
- If you find damage, stop using the device and contact OrthoVest for service or replacement parts.

STORAGE

- Hang the OrthoVest or lay it flat when not in use.
- Do not store it folded, crushed, or with the apron weight resting on it.
- Store in a clean, dry area away from direct heat.

LEAD APRON CARE

Follow your apron manufacturer's cleaning and inspection instructions. The OrthoVest does not replace required apron inspections.

CONTACT

Should you need any further information,
please do not hesitate to contact us.

123 Lincoln Blvd, Middlesex, NJ 08846, USA 

+1 866 575 1307 

service@phillips-safety.com 

www.phillips-safety.com 