

CLEANING & CARE GUIDELINES

INSTRUCTIONS OF USE

PHILLIPS-SAFETY.COM

COPYRIGHT ©2025. ALL RIGHTS RESERVED



PHILLIPS SAFETY
PRODUCTS INC.

USE

Phillips Safety Products must be used only for their designated purposes.

WARRANTY

Phillips Safety Radiation Apparel products are covered by a two-year warranty against manufacturing defects.

Attention: This does NOT include normal wear and overloading of use.

CLEANING

Use a mild detergent or cleaner to remove any dirt or debris.

Do NOT use abrasives.

DISINFECTING

It can be disinfected with most standard hospital-grade disinfectants containing QATs, IPA, or Hydrogen Peroxide.

CONTACT

Should you need any further information,
please do not hesitate to contact us.

123 Lincoln Blvd, Middlesex, NJ 08845, USA



+1 (866) 757 1307



service@phillips-safety.com



www.phillips-safety.com

